endovab

Endovascular Development AB

THE LIUNGMAN GUIDEWIRE FIXATOR

EN - Instruction for use	2
SV - Bruksanvisning	5
DE - Gebrauchsanweisung	8
FR - Mode d'utilisation	11
ES - Instruccions de uso	14

EN - Reference number
SV - Referensnummer
DE - Katalognummer
FR - Numéro de référence
ES - Número de referencia

EN - Do not use in MR fields
SV - Använd inte enheten i MR-fält
DE - Nicht in MR Feldern verwenden
FR - Ne pas exposer aux champs magnétiques
ES - No usar en campos de resonancia magnética

EN - Lot number
SV - Tillverkningsserie
DE - Chargennummer
FR - Numéro de lot
ES - Lote

EN - Sterilization with ethylene oxide
SV - Etylenoxidsteriliserad
DE - Sterilisation mit Ethylenoxide
FR - Stérilisation à l'oxyde d'éthylène
ES - Esterilización con óxido de Etileno

EN - Date of manufacture
SV - Tillverkningsdatum
DE - Herstellungsdatum
FR - Date de fabrication
ES - Fecha de manufactura

EN - Storage conditions
SV - Förvaringsförhållande
DE - Min. und max. Temperatur der Lagerung
FR - Conditions de stockage
ES - Condiciones de almacenamiento

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
THE LIUNGMAN GUIDEWIRE FIXATOR

DEVICE DESCRIPTION

The Liungman Guidewire Fixator (LGF) is a device for temporary fixation of a guidewire inside a blood vessel while maintaining the blood flow. The fixating element (A) is designed as a memory shaped Nitinol (Nickel Titanium alloy) mesh that, once deployed, produces a radial force against the inside of the vessel wall. The fixating element runs freely over the 0,035inch movable core guidewire (B) that is fitted with a stopper 30mm from the distal end, preventing the fixating element from traveling beyond the distal end. The guidewire comes in two lengths, 180cm and 280cm to support various applications. The fixating element is designed to expand radially and temporarily increase the radial force when moderate tension is applied to the guidewire.

The fixating element is delivered and retrieved through an introducer using the retrieval catheter (C) that is fitted with a spring shaped screw mechanism that connects to the thread cap on the fixating element. The torque device (D) can be attached to the retrieval catheter for better accuracy.

The Fixator is delivered inside an intubation sheath (E) to protect the fixating element when being inserted and retrieved through the introducer valve.

INTENDED USE

The Liungman Guidewire Fixator is intended to facilitate endovascular treatment methods where temporary fixation of a guidewire in a peripheral artery is indicated.

CONTRAINDICATIONS

Avoid positioning of the fixator in arteries with known risk for arterial dissection, or arteries with partial or complete arterial thrombosis.

WARNINGS

- Do not reuse
- The fixator is designed to withstand limited tension. Do not expose the fixated guidewire to excessive tension as this will increase the risk for arterial trauma.
- If incorrectly applied, guidewires can cause arterial perforations.
- The guidewires are fitted with distal cylindrical stoppers that increase the maximum guidewire diameter to 1,1mm (0,0433inch). Thus, test the combination of guidewire and medical device before use.
- Do not pull the guidewire upon retrieval as this may prevent fixator loading into the sheath.
- Do not over-rotate the retrieval catheter upon fixator retrieval.
- CAUTION: The device must not be used under magnetic resonance fields!

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in endovascular techniques. Before use, and if possible during the procedure, check the guidewire carefully for damage, kinks and other deformations. Damaged, misshapen, bent and otherwise deformed guidewires should not be used and should be discarded since there is a risk that the guidewire will be difficult to remove from the blood vessel.
- Confirmation of the position of the guidewire using fluoroscopy prevents accidental perforation of vessels during deployment of the device.
- Anticoagulant therapy sufficient to maintain an Activated Clotting Time of at least 250 seconds for the duration of the procedure is recommended.
- Regularly verify blood-flow over the fixator to reduce the risk for thrombus formations.
- Always use coaxial access in tortuous anatomies to allow smooth retrieval catheter rotation in the introducer sheath.

POTENTIAL ADVERSE EVENTS

- Hemorrhage
- Arterial dissection
- Arterial thrombosis
- Arterial embolism


INTRODUCER SHEATH SELECTION

Use Terumo Destination 7F, RSR06

INSTRUCTION FOR USE

Delivery

1. Measure the diameter (d) of the target artery and select appropriate Liungman Guidewire Fixator (LGF) size:
 - $d < 5\text{mm}$, LGF is not applicable
 - $5 \leq d \leq 7\text{mm}$, use 8mm LGF
 - $7 < d \leq 10\text{mm}$, use 11mm LGF
 - $10 < d \leq 13\text{mm}$, use 14mm LGF
 - $d > 13\text{mm}$, LGF is not applicable
2. Access the target artery and with a 7F introducer using conventional techniques.
3. Select the long or the short LGF guidewire for the intended application. Flush the guidewire dispenser completely with heparinized saline and then insert the guidewire through the introducer and position the guidewire in the target artery.
4. Remove the LGF retrieval catheter and the transparent intubation sheath containing the collapsed fixating element and flush the catheter completely with heparinized saline before removing the luer connection.
5. Make sure that the retrieval catheter is attached properly to the fixating element, then insert the retrieval catheter and fixating element over the guidewire into the introducer using the intubation sheath to access the valve. Note, there can be a slight resistance when inserting the guidewire through the fixator.
6. Push, and thus translate, the fixating element to the distal end of the introducer while using fluoroscopy to validate position in the target artery.
7. To deploy the fixating element, pull the introducer back over the retrieval catheter while maintaining the position of the retrieval catheter and fixating element at the desired position in the target artery.
8. Unscrew the retrieval catheter 2-3 rotations counter clockwise using the torque device and slowly retract and remove the retrieval catheter, while verifying under fluoroscopy that the fixator is detached. Use the transparent intubation sheath when exiting the hub to protect the retrieval mechanism. Also verify that the proximal end of the catheter doesn't catch and pull the movable core handle of the guidewire.

9. Optional: The guidewire has a removable stiff core that can be removed if a more pliable guidewire is

Retrieval

1. If the stiff core was removed, gently put it back into the guidewire without kinking it.
2. Insert the 7F introducer, with the dilator, over the LGF guidewire until positioned close to the fixating element.
3. Remove the dilator and insert the retrieval catheter over the guidewire into the introducer using the intubation sheath until the retrieval catheter reaches the distal end of the introducer.
4. Gently push the retrieval catheter towards the fixating element and screw the retrieval catheter clockwise 2-3 turns onto the fixating element using the torque device. Do not exceed the number of rotations as this can twist the fixator in the artery or in the sheath.

5. Slowly push the introducer while pulling the retrieval catheter to retrieve the fixator into the introducer.
6. Once collapsed inside the introducer, remove the retrieval catheter including the fixator completely out through the introducer using the intubation sheath when exiting the introducer valve.
7. When the position in the target artery is no longer needed, also remove the LGF guidewire.

Optional retrieval

If the normal retrieval procedure fails, the fixator can be retrieved using a conventional snare to catch the proximal end of the fixator and pull the fixator into a suitable larger introducer and retract the fixator by the snare.

STERILIZATION

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for single use only. Sterile if package is unopened and undamaged. Do not use the product if there is any breakage of the packaging. Upon removal from the package, inspect the product to ensure no damage has occurred.

STORAGE AND HANDLING

The products must be stored in a dry and dark place in the original box within a temperature range of 5–40°C (41–104°F).

REFERENCES

These instructions for use are based on experience from physicians and (or) published literature. Refer to your Endovab sales representative for information on available literature.

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