

Manual

PennyTalks (PT)


The personal digital voice labeler with speech supported user instructions and real-time backup.

Quick guide:

Put a self-adhesive label to whatever you would like to memorize.


Hold the sensor to the label and record your own voice memo in PT.

Playback your voice memo by holding the sensor of PT to the pre-recorded label.



Memory card is underneath the batteries:

PT is equipped with a removable memory card where all your memos are stored. This is an automatic backup of all your memos in case of any damage to the product.



3

How to set the language of the speech supported user instructions:

PennyTalks holds 12 pre-stored language modes for the speech supported user instructions.

A chime signal mode support function is available as well.

1. Press and hold the On/Off key for 3 seconds.
2. Touch to any of the languages available on the language setting chart with the sensor.
3. PT reports chosen language and which is now selected.

4

Switch on/switch off:

Insert the batteries. Press and hold the on/off button for 3 seconds. The speech support reports if PT is in an off or on mode.

If PT is left without any activity for more than 60 seconds, it will automatically switch off.

To record a label memo:

1. Touch the label with the sensor until PT reports "This is a new label".
2. Press and hold any of the + or - buttons until PT commands "Start recording".
3. Keep holding the button and record your memo.
4. When recording is done, release the button.
5. Your memo is now stored and PT reports "Recording is done"

Volume:

Press any of the + or - keys to raise or to lower the volume.

5

Playback:

Hold the sensor to the label. PT now reports your recorded memo.

Erase a recorded label memo:

1. Touch the label with the sensor.
2. PT starts to report your pre-recorded memo
3. Press and hold **BOTH** + and - keys until PT reports "Erased".

Actual label is now ready to be re-recorded.

Included in the box:

Alkaline batteries (not chargeable)
USB cable
Lanyard
Earphone
Self-adhesive tactile labels
Soft self-adhesive labels for clothing, safe to wash and to tumble dry.
Language setting chart

Technical specifications:

Batteries: Type AAA max 1.5 V
USB connection
Standard Memory size: up to 2700 memos
Optional memory size: up to 90000 memos