

# CAPSULAR TENSION RINGS

## Flexible Insertion, Stable Result

Highly flexible due to the compression molding technology. These PMMA Capsular Tension Rings are indicated for the stabilization of weakened, broken or missing zonules preventing IOL decentration after capsular shrinking.



### Main Features & Benefits

- ✓ **CIRCULAR EXPANSION AND STABILIZATION OF THE CAPSULAR BAG**
- ✓ **SAFE IOL CENTRATION IN EYES WITH ZONULAR DEHISCENCE**
- ✓ **PREVENTS IOL DECENTRATION AFTER CAPSULAR SHRINKING**
- ✓ **STABILIZED CONDITIONS DURING (PHACO EMULSIFICATION) SURGERY**
- ✓ **REDUCED RISK OF CAPSULAR FIBROSIS**

#### Capsular Tension Ring | Model 275 12/10




**Compression:** From 12 mm to 10 mm  
**Material:** PMMA | Flexible Ring  
**Available in:** Clear and brown PMMA  
**Overall Ø:** 12 mm

#### Capsular Tension Ring | Model 276 13/11



**Compression:** From 13 mm to 11 mm  
**Material:** PMMA | Flexible Ring  
**Available in:** Clear and brown PMMA  
**Overall Ø:** 13 mm



Recommended to be used in combination with the ON 385 EASY Control™ Micro Inserter