



# DOSI-FUSER®

	65ml	100ml	150ml	250ml	
7 days		<b>7days</b> 100ml 0.6ml/h 25915-100D7	<b>7days</b> 150ml 0.8ml/h 25915-150D7	<b>7days</b> 250ml 1.4ml/h 25915-250D7	Slow Flow with Capillary Element  
5 days	<b>5days</b> 65ml 0.5ml/h 25915-065D5	<b>5days</b> 100ml 0.8ml/h 25915-100D5	<b>5days</b> 150ml 1.2ml/h 25915-150D5	<b>5days</b> 250ml 2ml/h 25915-250D5	
2 days	<b>2days</b> 65ml 1.3ml/h 25915-065D2	<b>2days</b> 100ml 2.0ml/h 25915-100D2	<b>2days</b> 150ml 3.1ml/h 25915-150D2	<b>2days</b> 250ml 5.2ml/h 25915-250D2	
1 day	<b>1day</b> 65ml 2.7ml/h 25915-065D1	<b>1day</b> 100ml 4.1ml/h 25915-100D1	<b>1day</b> 150ml 6.2ml/h 25915-150D1	<b>1day</b> 250ml 10.4ml/h 25915-250D1	
12 hours	<b>12hours</b> 65ml 5.4ml/h 25915-065H12	<b>12hours</b> 100ml 8.3ml/h 25915-100H12	<b>12hours</b> 150ml 12ml/h 25915-150H12	<b>12hours</b> 250ml 20ml/h 25915-250H12	Fast Flow with Capillary Tube  
5 hours	<b>5hours</b> 65ml 13ml/h 25915-065H5	<b>5hours</b> 100ml 20ml/h 25915-100H5	<b>5hours</b> 150ml 30ml/h 25915-150H5	<b>5hours</b> 250ml 50ml/h 25915-250H5	
2 hours	<b>2hours</b> 65ml 32ml/h 25915-065H2	<b>2hours</b> 100ml 50ml/h 25915-100H2	<b>2hours</b> 150ml 75ml/h 25915-150H2	<b>2hours</b> 250ml 125ml/h 25915-250H2	
1 hours	<b>1hour</b> 65ml 65ml/h 25915-065H1	<b>1hour</b> 100ml 100ml/h 25915-100H1	<b>1hour</b> 150ml 150ml/h 25915-150H1	<b>1hour</b> 250ml 250ml/h 25915-250H1	
30 min.	<b>30 min</b> 65ml 130ml/h 25915-065M30	<b>30 min</b> 100ml 200ml/h 25915-100M30	<b>30 min</b> 150ml 300ml/h 25915-150M30		

Marknadsförs av Nordisk Medicin & Teknik AB

Box 3031, 434 03 Vallda

+46 300 32 30 83 [info@nmt-ab.se](mailto:info@nmt-ab.se)

[www.nmt-ab.se](http://www.nmt-ab.se)

